

Adenosine A1 Receptor Rabbit mAb


ABZ2579

Product Overview


Name	Adenosine A1 Receptor Rabbit mAb
Accession(Primary)	P30542
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:1000; ICC/IF-1:50-1:100; FC-1:100-1:1000; IP-1:20-1:100

Target information (P30542)

Gene Name	ADORA1
About Gene Symbol	



Western blot analysis of Adenosine A1 Receptor in K562, C6, 3T3, HeLa lysates using Adenosine A1 Receptor antibody.



Immunocytochemistry analysis of Adenosine A1 Receptor (green) in MCF-7 using Adenosine A1 Receptor antibody, and DAPI (blue)