

AdSS 2 Rabbit mAb


ABZ2580

Product Overview


Name	AdSS 2 Rabbit mAb
Accession(Primary)	P30520
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:10-1:100

Target information (P30520)

Gene Name	ADSS2
About Gene Symbol	ADSS2 (Adenylosuccinate Synthase 2) is a biosynthetic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of AdSS 2 in K562, rat Brain, C6, 3T3, HeLa lysates using AdSS 2 antibody.



Immunocytochemistry analysis of AdSS 2(green) in HeLa using AdSS 2 antibody, and DAPI(blue)