

B2R Rabbit mAb


ABZ2581

Product Overview

Name	B2R Rabbit mAb
Accession(Primary)	P30411
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 44 kDa; Observed MW: 80-100 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (P30411)

Gene Name	BDKRB2
About Gene Symbol	



Western blot analysis of BDKRB2 in A549, HL-60, U2OS lysates using BDKRB2 antibody.