

Aquaporin 1 Rabbit mAb


ABZ2587

Product Overview


Name	Aquaporin 1 Rabbit mAb
Accession(Primary)	P29972
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:200-1:2000

Target information (P29972)

Gene Name	AQP1
About Gene Symbol	AQP1 (Aquaporin 1 (Colton Blood Group)) is a formin actin nucleator involved in adaptive immunity. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Aquaporin 1 in C6 lysates using Aquaporin 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Aquaporin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.