

PML Protein Rabbit mAb


ABZ2588

Product Overview

Name	PML Protein Rabbit mAb
Accession(Primary)	P29590
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 98 kDa; Observed MW: 98,117 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:100-1:200; IP-1:20-1:50

Target information (P29590)

Gene Name	PML
About Gene Symbol	



Western blot analysis of PML Protein in A549
lysates using PML Protein antibody.