


VCAM1 Rabbit mAb

ABZ2589

Product Overview

Name	VCAM1 Rabbit mAb
Accession(Primary)	P29533
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 81 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P29533)

Gene Name Vcam1**About Gene Symbol**

Western blot analysis of VCAM1 in 3T3 lysates
using VCAM1 antibody.

