

Caspase 1 Rabbit mAb

ABZ2590


Product Overview

Name	Caspase 1 Rabbit mAb
Accession(Primary)	P29466
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45,20 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P29466)

Gene Name CASP1

About Gene Symbol



Western blot analysis of Caspase1 in HeLa lysates using Caspase 1 antibody.

