


CTGF Rabbit mAb

ABZ2593

Product Overview

Name	CTGF Rabbit mAb
Accession(Primary)	P29279
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 38 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P29279)

Gene Name CCN2**About Gene Symbol**

Western blot analysis of CTGF in C6, HeLa lysates using CTGF antibody.

