

TEF1/TEAD-1 Rabbit mAb


ABZ2596

Product Overview

Name	TEF1/TEAD-1 Rabbit mAb
Accession(Primary)	P28347
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:1000

Target information (P28347)

Gene Name	TEAD1
About Gene Symbol	TEAD1 (TEA Domain Transcription Factor 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TEF1/TEAD1 in C6, CHO-K1, HeLa lysates using Transcriptional Enhancer Factor TEF 1 antibody.