

POLD1 Rabbit mAb


ABZ2597

Product Overview



Name	POLD1 Rabbit mAb
Accession(Primary)	P28340
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 124 kDa; Observed MW: 124 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:100-1:200

Target information (P28340)

Gene Name	POLD1
About Gene Symbol	



Immunocytochemistry analysis of POLD1 (green) in K562 using POLD1 antibody, and DAPI (blue).



Western blot analysis of POLD1 in K562, C6, 3T3, HeLa lysates using POLD1 antibody.