

Choline Acetyltransferase Rabbit mAb


ABZ2598

Product Overview


Name	Choline Acetyltransferase Rabbit mAb
Accession(Primary)	P28329
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 83 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:100

Target information (P28329)


Gene Name	CHAT
About Gene Symbol	



Immunocytochemistry analysis of Choline Acetyltransferase (green) in SH-SY5Y using Choline Acetyltransferase antibody, and DAPI (blue).



Western blot analysis of Choline Acetyltransferase in mouse brain lysates using Choline Acetyltransferase antibody.



Western blot analysis of Choline Acetyltransferase
in rat Brain lysates using Choline Acetyltransferase
antibody