

LOX Rabbit mAb


ABZ2599

Product Overview

Name	LOX Rabbit mAb
Accession(Primary)	P28300
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:10-1:100


Target information (P28300)

Gene Name	LOX
About Gene Symbol	




Copyright© 2024 Zen-bioscience Limited.


Western blot analysis of LOX in HeLa, Jurkat lysates using LOX antibody.



Western blot analysis of LOX in (1) Jurkat lysates; (2) mouse brain lysates using LOX antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast using LOX antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of LOX in Jurkat using LOX antibody.