

PSMB8 Rabbit mAb


ABZ2601

Product Overview


Name	PSMB8 Rabbit mAb
Accession(Primary)	P28062
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 30 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:50-1:100; FC-1:50-1:100

Target information (P28062)

Gene Name	PSMB8
About Gene Symbol	PSMB8 (Proteasome 20S Subunit Beta 8) is a protease enzyme involved in ubiquitin-proteasome pathway and innate immune signaling. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate



Western blot analysis of LMP7 in 3T3, HeLa lysates using Proteasome beta 8 antibody.



Immunocytochemistry analysis of LMP7(green) in HeLa using LMP7 antibody, and DAPI(blue)