

PI3 Kinase p85 alpha Rabbit mAb

ABZ2602


Product Overview

Name	PI3 Kinase p85 alpha Rabbit mAb
Accession(Primary)	P27986
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P27986)

Gene Name PIK3R1

About Gene Symbol



Western blot analysis of PI 3 Kinase p85 alpha in HL-60, U2OS, C6 lysates using PI3 Kinase p85 alpha antibody.

