

PI3 Kinase p85 alpha Rabbit mAb


ABZ2603

Product Overview


Name	PI3 Kinase p85 alpha Rabbit mAb
Accession(Primary)	P27986
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:20000; IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (P27986)

Gene Name	PIK3R1
About Gene Symbol	



Western blot analysis of PI3 Kinase p85 alpha in Jurkat, Rat Brain, C6, CHO-K1 lysates using PI3 Kinase p85 alpha antibody.



Immunocytochemistry analysis of PI3 Kinase p85 alpha (green) in MCF-7 using PI3 Kinase p85 alpha antibody, and DAPI (blue)