

RPA70 Rabbit mAb


ABZ2606

Product Overview

Name	RPA70 Rabbit mAb
Accession(Primary)	P27694
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:10-1:50

Target information (P27694)

Gene Name	RPA1
About Gene Symbol	RPA1 (Replication Protein A1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of RPA70 in A549 lysates
using RPA70 antibody.