

HIF1 beta Rabbit mAb


ABZ2607

Product Overview


Name	HIF1 beta Rabbit mAb
Accession(Primary)	P27540
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; FC-1:200-1:500

Target information (P27540)


Gene Name	ARNT
About Gene Symbol	ARNT (Aryl Hydrocarbon Receptor Nuclear Translocator) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HIF1 beta in Raw264.7, 3T3 lysates using HIF1 beta antibody.



Western blot analysis of HIF1 beta in HeLa, CHO-K1, C6, Ramos lysates using HIF1 beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HIF1 beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.