


DPP4 Rabbit mAb

ABZ2608

Product Overview

Name	DPP4 Rabbit mAb
Accession(Primary)	P27487
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P27487)

Gene Name DPP4**About Gene Symbol**

Western blot analysis of CD26 in HeLa, A549, U251 lysates using CD26 antibody.

