


14-3-3 theta Rabbit mAb

ABZ2610


Product Overview

Name	14-3-3 theta Rabbit mAb
Accession(Primary)	P27348
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:200-1:1000

Target information (P27348)

Gene Name YWHAQ**About Gene Symbol**

Western blot analysis of 14-3-3 Theta in rat brain, rat lung lysates using 14-3-3 theta antibody.



Western blot analysis of 14-3-3 Theta/Tau in 3T3, HeLa lysates using 14-3-3 Theta/Tau antibody.