

MST1 Rabbit mAb


ABZ2612

Product Overview


Name	MST1 Rabbit mAb
Accession(Primary)	P26927
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 80 kDa; Observed MW: 96 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; FC-1:10-1:100

Target information (P26927)


Gene Name	MST1
About Gene Symbol	



Western blot analysis of MST1 in C2C12 lysates using MST1 antibody.



Western blot analysis of MST1 in HeLa, C6, rat Brain lysates using MST1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using MST1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.