

Phospholamban Rabbit mAb


ABZ2614

Product Overview

Name	Phospholamban Rabbit mAb
Accession(Primary)	P26678
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 6 kDa; Observed MW: 12,24 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; IP-1:20-1:50

Target information (P26678)

Gene Name	PLN
About Gene Symbol	



Western blot analysis of Phospho-lamban in mouse heart, rat heart lysates using Phospho-lamban antibody.