

PTBP1 Rabbit mAb


ABZ2615

Product Overview


Name	PTBP1 Rabbit mAb
Accession(Primary)	P26599
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:200-1:2000

Target information (P26599)


Gene Name	PTBP1
About Gene Symbol	PTBP1 (Polypyrimidine Tract Binding Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PTBP1 in 3T3 lysates using PTBP1 antibody.



Western blot analysis of PTBP1 in HeLa, CHO-K1, Ramos lysates using PTBP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PTBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.