

S100A4 Rabbit mAb


ABZ2617

Product Overview


Name	S100A4 Rabbit mAb
Accession(Primary)	P26447
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:100-1:500; IP-1:10-1:100

Target information (P26447)


Gene Name	S100A4
About Gene Symbol	



Western blot analysis of S100A4 in HeLa lysates using S100A4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using S100A4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of S1A4 (green) in
hela using S100A4 antibody, and DAPI(blue)