

S100A4 Rabbit mAb


ABZ2618

Product Overview

Name	S100A4 Rabbit mAb
Accession(Primary)	P26447
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1:20-1:50


Target information (P26447)

Gene Name	S100A4
About Gene Symbol	




Copyright© 2024 Zen-bioscience Limited.


Western blot analysis of S100A4 in A549 lysates using S100A4 antibody.



Western blot analysis of S100A4 in A375 lysates using S100A4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma using S100A4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of S100A4 in HeLa using S100A4 antibody.