

PAX6 Rabbit mAb


ABZ2619

Product Overview

Name	PAX6 Rabbit mAb
Accession(Primary)	P26367
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P26367)

Gene Name	PAX6
About Gene Symbol	PAX6 (Paired Box 6) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PAX6 in K562, HeLa
lysates using PAX6 antibody.