

Dnmt1 Rabbit mAb


ABZ2620

Product Overview


Name	Dnmt1 Rabbit mAb
Accession(Primary)	P26358
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 183 kDa; Observed MW: 183 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:100-1:1000; IP-1:20

Target information (P26358)

Gene Name	DNMT1
About Gene Symbol	DNMT1 (DNA Methyltransferase 1) is a methyltransferase involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Dnmt1 in Jurkat lysates using Dnmt1 antibody.



Immunocytochemistry analysis of Dnmt1(green) in HeLa using Dnmt1 antibody, and DAPI(blue)