

Moesin Rabbit mAb


ABZ2622

Product Overview


Name	Moesin Rabbit mAb
Accession(Primary)	P26038
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:20-1:50; ICC/IF-1:20-1:50; IP-1:20

Target information (P26038)


Gene Name	MSN
About Gene Symbol	MSN (Moesin) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Immunocytochemistry analysis of Moesin (green) in A549 using Moesin antibody, and DAPI (blue).



Western blot analysis of Moesin in C6, mouse heart lysates using Moesin antibody.



Western blot analysis of Moesin in A549, MCF-7, 293T, C2C12 lysates using Moesin antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Moesin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.