


LMO2 Rabbit mAb

ABZ2625

Product Overview

Name	LMO2 Rabbit mAb
Accession(Primary)	P25791
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (P25791)

Gene Name LMO2**About Gene Symbol**

Western blot analysis of LMO2 in Jurkat, HeLa, Raji lysates using LMO2 antibody.

