

PSMA4 Rabbit mAb


ABZ2626

Product Overview


Name	PSMA4 Rabbit mAb
Accession(Primary)	P25789
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:200-1:500; IP-1:100-1:200

Target information (P25789)

Gene Name	PSMA4
About Gene Symbol	



Western blot analysis of PSMA4 in K562, C6, 3T3, HeLa lysates using PSMA4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PSMA4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.