

ATP5A Rabbit mAb


ABZ2628

Product Overview

Name	ATP5A Rabbit mAb
Accession(Primary)	P25705
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:10000; IHC-1:50-1:2000; ICC/IF-1:100-1:200; FC-1:10-1:100


Target information (P25705)

Gene Name	ATP5F1A
About Gene Symbol	




Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of ATP5A in HepG2, HeLa lysates using ATP5A antibody.



Western blot analysis of ATP5A1 in HepG2 lysates; using ATP5F1A antibody.



Western blot analysis of ATP5F1A in HeLa □Rat brain□Mouse brain lysates using ATP5F1A antibody.