

BRD2 Rabbit mAb


ABZ2631

Product Overview

Name	BRD2 Rabbit mAb
Accession(Primary)	P25440
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 88 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC:1:20-1:50; ICC/IF-1:20-1:50; FC-1:20-1:50

Target information (P25440)

Gene Name	BRD2
About Gene Symbol	



Western blot analysis of BRD2 in rat Brain, Jurkat lysates using BRD2 antibody.