

Myelin Protein Zero Rabbit mAb


ABZ2632

Product Overview

Name	Myelin Protein Zero Rabbit mAb
Accession(Primary)	P25189
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; IP-1:20

Target information (P25189)

Gene Name	MPZ
About Gene Symbol	



Western blot analysis of Myelin Protein Zero in C6
lysates using Myelin Protein Zero antibody.