


RDC1/CXCR-7 Rabbit mAb

ABZ2633

Product Overview

Name	RDC1/CXCR-7 Rabbit mAb
Accession(Primary)	P25106
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P25106)

Gene Name ACKR3**About Gene Symbol**

Western blot analysis of GPCR RDC1/CXCR7 in K562, C6, 3T3, HeLa lysates using RDC1 antibody.

