


Neurokinin 1 Receptor Rabbit mAb

ABZ2634


Product Overview

Name	Neurokinin 1 Receptor Rabbit mAb
Accession(Primary)	P25103
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:50-1:100

Target information (P25103)

Gene Name TACR1**About Gene Symbol**

Western blot analysis of Neurokinin 1 Receptor in 3T3 lysates using Neurokinin 1 Receptor antibody.



Western blot analysis of Neurokinin 1 Receptor in CHO-K1, C6, Jurkat lysates using Neurokinin 1 Receptor antibody.