

APC Rabbit mAb


ABZ2636

Product Overview



Name	APC Rabbit mAb
Accession(Primary)	P25054
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 312 kDa; Observed MW: 160 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; IHC-1:50-1:100

Target information (P25054)

Gene Name	APC
About Gene Symbol	



Immunocytochemistry analysis of APC (green) in LNCap using APC antibody, and DAPI (blue).



Western blot analysis of APC in HeLa, A549, HL-60, U2OS lysates using APC antibody.