

Galectin 3 Rabbit mAb


ABZ2639

Product Overview


Name	Galectin 3 Rabbit mAb
Accession(Primary)	P17931
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:500-1:1000; FC-1:20-1:50

Target information (P17931)


Gene Name	LGALS3
About Gene Symbol	LGALS3 (Galectin 3) is an extracellular matrix laminin protein involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Galectin 3 in rat lung lysates using Galectin 3 antibody.



Western blot analysis of Galectin 3 in 3T3, HeLa lysates using Galectin 3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Galectin 3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.