

Galectin 3 Rabbit mAb


ABZ2640

Product Overview


Name	Galectin 3 Rabbit mAb
Accession(Primary)	P17931
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:20-1:100

Target information (P17931)


Gene Name	LGALS3
About Gene Symbol	LGALS3 (Galectin 3) is an extracellular matrix laminin protein involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Galectin 3 (green) in U87-MG using Galectin 3 antibody, and DAPI(blue).



Western blot analysis of Galectin 3 in rat colon, rat lung lysates using Galectin 3 antibody.



Western blot analysis of Galectin 3 in HeLa, COS7
lysates using Galectin 3 antibody