

Galectin 3 Rabbit mAb


ABZ2641

Product Overview


Name	Galectin 3 Rabbit mAb
Accession(Primary)	P17931
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:1000; HC-P:1:200-1:2000; ICC/IF-1:200-1:500; FC-1:20-1:50; IP-1:20-1:50; IHC-1:50-1:100

Target information (P17931)


Gene Name	LGALS3
About Gene Symbol	LGALS3 (Galectin 3) is an extracellular matrix laminin protein involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Galectin 3 (green) in U87-MG using Galectin 3 antibody, and DAPI(blue).



Western blot analysis of Galectin 3 in HeLa lysates using Galectin 3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using Galectin 3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.