

Hsp70 Rabbit mAb


ABZ2642

Product Overview


Name	Hsp70 Rabbit mAb
Accession(Primary)	P17879
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:20-1:50; FC-1:100-1:200

Target information (P17879)

Gene Name	Hspa1B
About Gene Symbol	



Western blot analysis of Hsp70 in K562, rat Brain, 3T3, HeLa lysates using Hsp70 1B antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Hsp70 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.