


XBP1 Rabbit mAb

ABZ2643


Product Overview

Name	XBP1 Rabbit mAb
Accession(Primary)	P17861
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:50-1:100; IP-1:20-1:50

Target information (P17861)

Gene Name XBP1**About Gene Symbol**

Western blot analysis of XBP1 in Raw264.7
lysates using XBP1 antibody.



Western blot analysis of XBP1 in HL-60, U2OS
lysates using XBP1 antibody.