

PFKL Rabbit mAb


ABZ2644

Product Overview

Name	PFKL Rabbit mAb
Accession(Primary)	P17858
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:10-1:100

Target information (P17858)

Gene Name	PFKL
About Gene Symbol	



Western blot analysis of PFKL in C6, 3T3 lysates using PFKL antibody.