

DDX5 Rabbit mAb


ABZ2645

Product Overview


Name	DDX5 Rabbit mAb
Accession(Primary)	P17844
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 69 kDa; Observed MW: 69 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1:10-1:100

Target information (P17844)


Gene Name	DDX5
About Gene Symbol	DDX5 (DEAD-Box Helicase 5) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of DDX5 in C6, 3T3, HeLa lysates using DDX5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using DDX5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of DDX5 (green) in Hela using DDX5 antibody, and DAPI (blue)