

HOXB4 Rabbit mAb

ABZ2649

Product Overview


Name	HOXB4 Rabbit mAb
Accession(Primary)	P17483
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 28 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (P17483)

Gene Name	HOXB4
About Gene Symbol	HOXB4 (Homeobox B4) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of HOXB4 (green) in K562 using HOXB4 antibody, and DAPI (blue).



Western blot analysis of HOXB4 in K562 lysates using HOXB4 antibody.