

Spi1 Rabbit mAb


ABZ2650

Product Overview


Name	Spi1 Rabbit mAb
Accession(Primary)	P17433
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:20-1:50; IP-1:20-1:50

Target information (P17433)


Gene Name	Spi1
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse liver using SPI1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PU.1/Spi1 in CHO-K1, C6 lysates using PU.1/Spi1 antibody.



Immunocytochemistry analysis of PU.1/Spi1
(green) in Raw264.7 using PU.1/Spi1 antibody, and
DAPI(blue)