

Integrin alpha 2 Rabbit mAb


ABZ2651

Product Overview

Name	Integrin alpha 2 Rabbit mAb
Accession(Primary)	P17301
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 129 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200

Target information (P17301)

Gene Name	ITGA2
About Gene Symbol	ITGA2 (Integrin Subunit Alpha 2) is an integrin cell adhesion receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and iden



Western blot analysis of Integrin alpha 2 in HeLa lysates using Integrin alpha 2 antibody.