

JunB Rabbit mAb


ABZ2652

Product Overview

Name	JunB Rabbit mAb
Accession(Primary)	P17275
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 42,43 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100-1:200; IP-1:10-1:100

Target information (P17275)

Gene Name	JUNB
About Gene Symbol	JUNB (JunB Proto-Oncogene, AP-1 Transcription Factor Subunit) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of JunB in K562, C6, 3T3, HeLa lysates using JunB antibody.