

Interferon alpha/beta Receptor 1 Rabbit mAb


ABZ2653

Product Overview

Name	Interferon alpha/beta Receptor 1 Rabbit mAb
Accession(Primary)	P17181
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 64 kDa; Observed MW: 110-130 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P17181)

Gene Name	IFNAR1
About Gene Symbol	



Western blot analysis of Interferon alpha/beta receptor 1 in K562, HeLa lysates using Interferon alpha/beta Receptor 1 antibody.