

## HMGA1 Rabbit mAb


ABZ2655

### Product Overview


|                             |                                                                                           |
|-----------------------------|-------------------------------------------------------------------------------------------|
| <b>Name</b>                 | HMGA1 Rabbit mAb                                                                          |
| <b>Accession(Primary)</b>   | P17096                                                                                    |
| <b>Host Species</b>         | Rabbit                                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                                |
| <b>Isotype</b>              | IgG                                                                                       |
| <b>Species Reactivity</b>   | Human, Mouse                                                                              |
| <b>Calculated MW</b>        | Calculated MW: 12 kDa; Observed MW: 18 kDa                                                |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:1000-1:5000;<br>ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:20-1:50 |

### Target information (P17096)


|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | HMGA1 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of HMGA1 in HeLa, C6, MCF-7 lysates using HMGA1 antibody.



Immunocytochemistry analysis of HMGA1(green) in HeLa using HMGA1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using HMGA1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.