

ATF2 Rabbit mAb


ABZ2656

Product Overview


Name	ATF2 Rabbit mAb
Accession(Primary)	P16951
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; FC-1:20-1:50; IP-1:20-1:50

Target information (P16951)

Gene Name	Atf2
About Gene Symbol	



Western blot analysis of ATF2 in 3T3, PC-12 lysates using ATF2 antibody.



Western blot analysis of ATF2 in Jurkat lysates using ATF2 antibody.