

SERCA2 ATPase Rabbit mAb


ABZ2657

Product Overview


Name	SERCA2 ATPase Rabbit mAb
Accession(Primary)	P16615
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 115 kDa; Observed MW: 115 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:100-1:1000

Target information (P16615)


Gene Name	ATP2A2
About Gene Symbol	ATP2A2 (ATPase Sarcoplasmic/Endoplasmic Reticulum Ca ²⁺ Transporting 2) is a membrane transporter involved in membrane transport. It is an important research target in cardiovascular disease and is widely studied to elucidate cellular mechanisms and identi




Western blot analysis of SERCA2 ATPase in rat brain lysates using SERCA2 ATPase antibody.



Western blot analysis of SERCA2 ATPase in K562, HeLa lysates using SERCA2 ATPase antibody.



Immunocytochemistry analysis of SERCA2 ATPase(green) in Hela using SERCA2 ATPase antibody,and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using SERCA2 ATPase antibody.High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.